

Jerry Ford Company, LLC

Manual Wheelchairs with Fall Intervention & Safety Systems

Visit Us On The Web At:
www.JerryFordCompany.com



Seeing our innovative and effective brake systems in action will make you a believer!

Visit Us On-line at:
www.JerryFordCompany.com to view our "must see" video demonstration or to order our award-winning automatic brake systems!



Order On-line at:
www.JerryFordCompany.com

Or Email us at:
info@JerryFordCompany.com



Jerry Ford Company, LLC
29309 Hwy 74
St. Charles, MN 55972

P: 507-932-8866
F: 507-932-8833
TF: 866-800-6049

**"DOING ALL YOU
CAN IS PEACE OF MIND!"**



**PREVENT FALLS &
INJURIES WITH THE**

Safe-T-Chair™

Featuring
Jerry-Lock Brake Systems
Incline Brake Systems
Attendant Brake Systems



**Featuring Quality
Everest Jennings Wheelchairs**

Introducing the Award-Winning *Safe-T-Chair™*

Get one of our award-winning Jerry-Lock brake systems pre-mounted and ready to use on quality fixed arm wheelchairs *starting at just \$450.00*


JERRY-LOCK BRAKE SYSTEM

Automatically applies wheel brakes whenever the occupant is transferring to or from the wheelchair ensuring that the wheelchair is in a safe and stable “locked” condition. Manual Brake lock-offs allow the wheelchair to be easily brought along when occupant is ambulating.



INCLINE BRAKE SYSTEM

Helps attendant control the wheelchair on inclines, and prevents attendant strain and injury by providing the attendant with a hand lever friction brake.



ATTENDANT BRAKE SYSTEM

For Hospital escort or for wheelchair occupants whose safety may be enhanced by limiting brake control to only the attendant, Jerry Ford Company offers the Attendant Brake System.



- Truly user-friendly Automatic Wheelchair Brake System.
- Palm Release for Easy Movement of Unoccupied Wheelchair
- Lock-off Mechanism for Ambulating Patients
- Lock-off Auto Reset for Added Safety
- No Under-seat Pressure Points
- Does Not Interfere with Normal Folding or Operation

“Our automatic wheelchair brakes prevent falls and injuries, free up staff, and improve quality of life... an inexpensive investment that saves money, time and lives!”

SIZE OPTIONS

The Jerry-Lock Brake System and/or Incline Brake System can be ordered premounted on full length, fixed arm wheelchairs, in widths of 16”, 18”, or 20”.*

The Attendant Brake System and/or Incline Brake System can be ordered premounted on the full or desk length, removable arm wheelchair, in widths of 16”, 18”, or 20”.*

We also offer these systems on heavy duty wheelchairs with weight capacities of 450 lbs that come in 22”, 24”, or 26”

*Standard 16-18” chair weight capacity: 250 lbs, 20”: 300 lbs

www.JerryFordCompany.com Offers

Online Ordering, Video Demos, In-the-News, and Customer Testimonials

Jerry Ford Company, LLC

29309 Hwy 74

St. Charles, MN 55972

Call between 9am-3pm CST to speak with a customer representative!

TF: 866-800-6049

F: 507-932-8833

E: info@JerryFordCompany.com